

The Evolution of the Hologic SecurView® Workstation

At Hologic, we strive to bring you the best technology and workflow to assist you in providing optimum patient care. This review of the SecurView Standalone and Client workstation hardware shows the technology advancements we have made in recent years.

	Operating System (shipped configuration)	RAM	Disk Drive Redundancy	Encrypted Hard Drives	Study Storage Area	Meets minimum specifications for 3D Mammography™
Precision T5820 (5/2020)	Windows 10	32GB	Yes	Yes	7TB	Yes
Precision T3620 (12/2018)	Windows 7 Windows 10	32GB	Yes	Yes	7TB	Yes
Precision T7910 (6/2015)	Windows 7	16GB	Yes	Yes	3.5TB/7TB	Yes
Precision T7600/T7610 (12/2013)	Windows 7	16GB	Yes	No	3.5TB	Yes
Precision T7500 (6/2010)	Window XP Windows 7	6GB	Yes	No	1.6TB	Yes (Requires Windows 7)
Precision T7400 (7/2008)	Windows XP	4GB	No	No	976GB	No


If you would like to upgrade your current system in order to take advantage of the higher performance and increased capacity of the latest model, please contact your Hologic representative.

The Evolution of the Hologic SecurView® Workstation

At Hologic, we strive to bring you the best technology and workflow to assist you in providing optimum patient care. This review of the SecurView Standalone and Client workstation hardware shows the technology advancements we have made in recent years.

	Operating System (shipped configuration)	Processor	RAM	Study Storage Area	Disk Drive Redundancy	Encrypted Hard Drive	Dual Network Card & Dual Power Supply	Meets minimum specifications for 3D Mammography™
PowerEdge T640 (Tower) (1/2021)	Windows Server 2016	Intel Xeon Gold 5120T 2.2 GHz	64GB	40TB	Yes	Yes	Yes	Yes
PowerEdge T630 High Performance (Tower) (2/2018)	Windows Server 2016	Intel Xeon E5-2620 v4	64GB	24TB	Yes	Yes	Yes	Yes
PowerEdge T630 (Tower) PowerEdge R730 (Rackmount) (6/2015)	Windows Server 2008 and Windows Server 2016	Intel Xeon E5-2620 v3	16GB	5TB	Yes	Yes	Yes	Yes
PowerEdge T620 (Tower) PowerEdge R720 (Rackmount) (5/2012)	Windows Server 2008	Intel Xeon E5-2620	16GB	5TB	Yes	No	Yes	Yes
PowerEdge T610 (Tower) (12/2010)	Windows Server 2003 and Windows Server 2008	Intel Xeon E5620	6GB	3.6TB	Yes	No	3.5TB	Yes

If you would like to upgrade your current system in order to take advantage of the higher performance and increased capacity of the latest model, please contact your Hologic representative.

 2797 EC REP Hologic BV, Da Vincilaan 5, 1930 Zaventem, Belgium

Breast & Skeletal Solutions | www.3dimensionsmammography.eu | +32.2.711.4680